

# Orion Back




The Orion Back™ was designed for users who require mild to moderate postural support. Extras, such as the positioning options, the ease of attachment and the unique method of adjusting for a wide range of angles, makes the Orion Back appealing for a wide variety of users.


Prism Series Backs		STYLE: Orion Back	
Stock Sizes – Same Day Shipment			
HCPCS CODE: E2611			
Part Number	Size W x H	Part Number	Size W x H
OB1616	16" x 16"	OB1820	18" x 20"
OB1620	16" x 20"	OB2016	20" x 16"
OB1816	18" x 16"	OB2020	20" x 20"
Non-Stock Sizes – Shipment within 3 business days			
HCPCS CODE: E2611			
Part Number	Size W x H	Part Number	Size W x H
OB1414	14" x 14"	OB1818	18" x 18"
OB1618	16" x 18"	OB1916	19" x 16"
OB1716	17" x 16"	OB1920	19" x 20"
OB1720	17" x 20"		

## Product Highlights:

- ✓ Lightweight, sturdy 5052 grade aluminum shell
- ✓ A full 5" range of forward and rear adjustable
- ✓ Fully Adjustable for height, depth and angle of recline
- ✓ Attractive cover design
- ✓ Adjustable Angle Hardware
- ✓ Replaceable foam
- ✓ **Product Weight:** 4.75 lbs, (18"W x 16"H)
- User Weight:** 250 lbs
- Warranty:** 3 years limited on back, 90 days on cover, and 1 year on hardware



Lateral Dimensions		- Orion Back
Heights 11-26"		
Chair Width	DIM B (")	DIM A (")
14"	11.85	5
16"	13.85	5
18"	15.85	5
20"	17.85	5
22"	19.85	5



Cover and Foam Layers

✓ Hardware fits 3/4", 7/8" and 1" tubing and can be installed within a few minutes with the right tools