

Serene Champion



HCPCS: E2603/04 Skin Protection Cushion

THE CHAMPION CUSHION was designed for users that require mild postural support and offers gentle contours helping to promote correct seating posture. The cushion has moderate leg troughs and abductor that provides pelvic stability and leg positioning. A cushion made from medium density, high resilient, polyurethane foam with a 1" visco insert which provide additional support to the IT area.

PRODUCT HIGHLIGHTS:

- ✓ Standard water proof cover
- ✓ radiused front edge for foot propelling
- ✓ Light weight cushion
- ✓ Non slip bottom on washable, incontinent cover
- ✓ Promotes mild postural support
- ✓ 4 way stretch, water resistant inner liner
- ✓ Medium density, high resilient, polyurethane foam
- ✓ High density Visco Foam Insert


Product Weight: 3 Lbs / 1.36 Kg (20"W x 20"L)




Non slip Bottom on washable, incontinent cover



Polyurethane foam, anatomically shaped



Pressure Reading Scale in mmHG



User Weight 250 Lbs / 113 Kg
 Duration 1 hr
 Average Pressure 32 mmHg

SERENE CUSHION		STYLE: CHAMPION	
Stock Sizes – Same Day Shipment			
Part #	Size W x L	Part #	Size W x L
CH1616	16" x 16"	CH1818	18" x 18"
CH1618	16" x 18"	CH2016	20" x 16"
CH1816	18" x 16"	CH2018	20" x 18"

SERENE CUSHION		STYLE: CHAMPION	
Non-Stock Sizes – Shipment within 3 business days			
Part #	Size W x L	Part #	Size W x L
CH1620	16" x 20"	CH2216	22" x 16"
CH1820	18" x 20"	CH2218	22" x 18"
CH2020	20" x 20"	CH2220	22" x 20"
CH2416	24" x 16"	CH2418	24" x 18"
CH2420	24" x 20"		