

Account Information		Order Date: _____	
Account #: _____		P.O.# _____	
Business Name: _____		Phone: _____	
Ship to Address: _____		City: _____	
Fax: _____	Province: _____	Zip Code: _____	
Contact: _____			

**The HCPCS codes provided are based on code verification by the PDAC or our interpretation of the code definitions.
Proper use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.
Prices are subject to change without notice**

① FRAME TYPE				Standard Chair Includes:
<input checked="" type="checkbox"/>	Superchair for kids	104 E1233	\$2,765	
<input type="checkbox"/>	10"-13" Frame (crosst tubing)	104CT1	No Charge	
<input type="checkbox"/>	13"-16" Frame (crosst tubing)	104CT2	No Charge	
② FRAME FINISH				Rigid Tilting Frame w/Swingaway Front End Round Heavyduty Crossbrace Tubing Seat Width 10"-16" Seat Depth 11"-16" Seat to Floor Height Front 17.50"-20" Seat to Floor Height Rear 15.75"-18.5" Rear Wheel Options 12", 20" or 22") Handrim Option (Stainless Steel or Plastic Coated) Caster Options 4", 5", 7" or 8" Adjustable Height Full Length Arm (vinyl pad) T-Style Angle Adjustable 10 deg Bend Cane w/Push bar - High 70 Degree Swingaway Front Riggings Composite Flip-Up Footplate Wheel Locks (Push or Pull to Lock) Side Storage Pocket Padded Headrest (Kids or Small) Choice of Seven Colours (Black, Blue, Burgundy, Orange, Yellow, Red or Silver)
<input type="checkbox"/>	Black	104C01	No Charge	
J <input type="checkbox"/>	Blue	104C03	No Charge	
J <input type="checkbox"/>	Burgundy	104C11	No Charge	
J <input type="checkbox"/>	Orange	104C05	No Charge	
J <input type="checkbox"/>	Yellow	104C08	No Charge	
J <input type="checkbox"/>	Red	104C07	No Charge	
J <input type="checkbox"/>	Silver	104C02	No Charge	
③ TRANSPORT TIE DOWN				
<input type="checkbox"/>	Transport Tie Down System	104MA71	No Charge	
<input type="checkbox"/>	Omit Transport Tie Down System	104MA72	-\$150	

Transport Tie Down System includes four factory installed wheelchair transport brackets. Use these transport brackets only to secure an unoccupied wheelchair during transport. The Department of Transportation has not approved any tie-down systems for transportation of a user while sitted in a wheelchair, in a moving vehicle of any type. It is recommended that wheelchair users be transported in appropriate seating in moving vehicles and the vehicles restraints made available by the auto industry be used. Future Mobility Healthcare does not recommend any wheelchair transportation systems.



FOOTNOTES:
A. Seat angle is set at 3 degrees. See service manual for details
J. Allow a few extra days for powder coating time for custom colors

③ SEAT WIDTH			
<input type="checkbox"/>	10" Seat Width	104SW10	No Charge
<input type="checkbox"/>	11" Seat Width	104SW11	No Charge
<input type="checkbox"/>	12" Seat Width	104SW12	No Charge
<input type="checkbox"/>	13" Seat Width	104SW13	No Charge
<input type="checkbox"/>	14" Seat Width	104SW14	No Charge
<input type="checkbox"/>	15" Seat Width	104SW15	No Charge
<input type="checkbox"/>	16" Seat Width	104SW16	No Charge
④ SEAT DEPTH			
<input type="checkbox"/>	11" Seat Depth	104SD11	No Charge
<input type="checkbox"/>	12" Seat Depth	104SD12	No Charge
<input type="checkbox"/>	13" Seat Depth	104SD13	No Charge
<input type="checkbox"/>	14" Seat Depth	104SD14	No Charge
<input type="checkbox"/>	15" Seat Depth	104SD15	No Charge
<input type="checkbox"/>	16" Seat Depth	104SD16	No Charge

SUPERCHAIR FOR KIDS FRONT STF HEIGHT AVAILABILITY									
REAR WHEEL		4" CASTER							
12"									
20"	17.5"	18"							
22"									
REAR WHEEL		5" CASTER							
12"			18.5"						
20"									
22"			18.5"						
REAR WHEEL		7" CASTER							
12"				19"	19.5"				
20"									
22"				19"	19.5"				
REAR WHEEL		8" CASTER							
12"							20"		
22"							19.75"		

5 SEAT-TO-FLOOR HEIGHT - front - 20", 22" Rear Wheels

- 17.5" Seat to Floor (20" Wheels and 4" Casters) 104STF62 No Charge
- 18" Seat to Floor (20" Wheels and 4" Casters) 104STF63 No Charge
- 18.5" Seat to Floor (22" Wheels and 5" Casters) 104STF64 No Charge
- 19" Seat to Floor (22" Wheels and 7" Casters) 104STF60 No Charge
- 19.5" Seat to Floor (22" Wheels and 7" Casters) 104STF61 No Charge
- 19.75" Seat to Floor (22" Wheels and 8" Casters) 104STF65 No Charge
- 18.5" Seat to Floor (12" Wheels and 5" Casters) 104STF66 No Charge
- 19" Seat to Floor (12" Wheels and 7" Casters) 104STF67 No Charge
- 19.5" Seat to Floor (12" Wheels and 7" Casters) 104STF68 No Charge
- 20" Seat to Floor (12" Wheels and 8" Casters) 104STF69 No Charge

6 REAR WHEEL RIM

- Mag Wheels----- 104RW09 No Charge
- Spoke Wheels----- 104RW10 \$75

8 REAR WHEEL TIRE

- Airless (Solid Rubber)----- 104RT1 Standard
- Pneumatic----- 104RT2 E2211 \$90

7 HANDRIM OPTIONS

- Chrome----- 104HR1 Standard
- Plastic Coated----- 104HR2 K0059 \$160
- Omit Handrims----- 104HR3 less net \$15

FOOTNOTES:
I. Available in Desk length or Full Length Vinyl

8 ARM TYPE

- T-Style Adjustable Height Pin Release - Youth (vinyl)----- 104AR05 E0973 Standard
- T-Style Adjustable Height Pin Release - Youth (gel)----- 104AR14 E0973 \$205
- T-Style Fixed Height ----- 104AR02 E0973 No Charge
- Cantilever Arm Adj. Angle ----- 104AR11 E0973 \$225
- Cantilever Arm Adj. Angle with Guard ----- 104AR12 E0973 \$350

9 ARM PAD

- Desk length Vinyl 10" (pair)----- 104AP1 No Charge
- Full length Vinyl 14" (pair)----- 104AP2 No Charge
- Contour Gel 17" (pair)----- 104AP3 Standard

10 BACK TYPE

- Angle Adj 10 deg Bend Cane w/Push bar (low)----- 104SP3 No Charge
- Angle Adj 10 deg Bend Cane w/Push bar (high)----- 104SP4 Standard
- Angle Adj 10 deg Bend Cane w/o Push bar (low)----- 104PC3 No Charge
- Angle Adj 10 deg Bend Cane w/o Push bar (high)----- 104PC4 No Charge

11 BACK OPTIONS

- Super Chair Kids Back Stock ----- 104BK75 E2620 \$695
- Super Chair Kids Back Non Stock ----- 104BK76 E2620 \$695
- Super Chair Kids Back Custom ----- 104BK77 E2620 \$745
- Slip on Padded Back Upholstery----- 104BK48 \$150
- Prism Basic Back on chair ----- 104BK04 E2620 \$615
- Prism Ultra Back on chair ----- 104BK09 E2620 \$615
- Prism Air Back on chair ----- 104BK14 E2615 \$655
- Prism TrueFitt Back on chair ----- 104BK19 E2615 \$625
- Prism Max Ultra Back on chair ----- 104BK24 E2620 \$645
- Prism Max Air Back on chair ----- 104BK29 E2620 \$725
- Prism Max TrueFitt Back on chair ----- 104BK34 E2620 \$675
- Omit Back Upholstery----- 104BK49 Standard
- Specify Back Size W____x H____
- Specify Max Back Colour _____

12 BACK PINS (Prism Backs only)

Longer Pins for Back Hardware 1" Growth -----	104BP1	\$50
Longer Pins for Back Hardware 2" Growth -----	104BP2	\$50

13 SEAT OPTIONS

J	<input type="checkbox"/> Expandable Seat Pan-----	104SE1	E2231	\$150
	<input type="checkbox"/> Super Chair Kids Cushion Stock-----	104SC44	E2607	\$365
	<input type="checkbox"/> Super Chair Kids Cushion Non Stock-----	104SC45	E2607	\$365
	<input type="checkbox"/> Super Chair Kids Cushion Custom-----	104SC46	E2607	\$450
	<input type="checkbox"/> Prism Ideal Cushion -----	104SC03	E2601	\$100
	<input type="checkbox"/> Prism Elite Cushion -----	104SC06	E2603	\$225
	<input type="checkbox"/> Prism Supreme Cushion -----	104SC09	E2607	\$365
	<input type="checkbox"/> Prism Supreme Air Cushion -----	104SC15	E2624	\$475
	<input type="checkbox"/> Prism Supreme II Plus Cushion -----	104SC12	E2607	\$385
	<input type="checkbox"/> Prism Drop Base on chair FM -----	104SC36	E2231	\$100
	<input type="checkbox"/> Prism Drop Base on chair 1" -----	104SC37	E2231	\$100
	<input type="checkbox"/> Prism Drop Base on chair 2" -----	104SC38	E2231	\$100
	<input type="checkbox"/> Prism Adjustable Drop Base on chair -----	104SC42	E2231	\$150
B	<input type="checkbox"/> Omit Cushion-----	104SC20		less net \$10

FOOTNOTES:

- B. Super Chair Kids Aluminum Footplate available for chairs 10"W to 16" W
- C. Composite Flip-Up Footplates available for chairs 14" to 16" W
- D. Aluminum Flip-Up Footplates available for chairs 12" to 16" W
- E. Aluminum Angle Adj. Flip-Up Footplates available for chairs 15"W to 16"W
- J. Expanable seat pan available in depths 11-14 and 13-16

14 FRONT RIGGING

C	<input type="checkbox"/> 70 Degree Swingaway Upgraded-----	104FR04		Standard
D	<input type="checkbox"/> 80 Degree Swingaway Upgraded-----	104FR19		\$97
E	<input type="checkbox"/> 90 Degree Swingaway Upgraded-----	104FR05		\$97
	<input type="checkbox"/> Elevating Legrest with Swingaway Upgraded-----	104FR10	E0990	\$250
	<input type="checkbox"/> Omit Front Rigging-----	104FR11		less net \$20

15 FOOTPLATE TYPE

	<input type="checkbox"/> Super Chair Kids Aluminum Flip-Up Footplate Small-----	104FP13		No Charge
	<input type="checkbox"/> Super Chair Kids Aluminum Flip-Up Footplate Medium-----	104FP14		No Charge
	<input type="checkbox"/> Super Chair Kids Aluminum Flip-Up Footplate Large-----	104FP15		No Charge
H	<input type="checkbox"/> Composite Flip-Up Footplate Small-----	104FP01		Standard
	<input type="checkbox"/> Aluminum Flip-Up Footplate Medium-----	104FP05		\$48
	<input type="checkbox"/> Aluminum Adjustable Angle Flip-Up Footplate Small-----	104FP10	K0040	\$159

16 FOOTREST/LEGREST ACCESSORIES

	<input type="checkbox"/> Padded Calf Strap-----	104FA03	K0038	\$35
	<input type="checkbox"/> Foam Pad on Footplate (pair)-----	104FA07		\$70
	<input type="checkbox"/> Heel Loops-----	104FA01	E0951	\$40
	<input type="checkbox"/> Omit Heel Loops-----	104FA02		less net \$7

17 WHEEL LOCKS

	<input type="checkbox"/> Push-to-Lock Wheel Lock-----	104WL1		Standard
	<input type="checkbox"/> Pull-to-Lock Wheel Lock-----	104WL2		No Charge
	<input type="checkbox"/> Wheel Lock Extension-----	104WL4	E0961	\$75

18 HEADREST HARDWARE OPTIONS

	<input type="checkbox"/> Fixed Headrest Hardware -----	104HE01	E0955	Standard
	<input type="checkbox"/> Single Multi-Axial Headrest Hardware -----	104HE03	E0955	\$217
F.G	<input type="checkbox"/> Dual Multi-Axial Headrest hardware -----	104HE05	E0955	\$219
	<input type="checkbox"/> Fixed Headrest Hardware Ass'y (+ 4")-----	104HE02	E0955	\$20
	<input type="checkbox"/> Single Multi-Axial Headrest Hardware (+4") -----	104HE04	E0955	\$222
	<input type="checkbox"/> Dual Multi-Axial Headrest hardware (+ 4") -----	104HE06	E0955	\$225
	<input type="checkbox"/> Links Multi Adjustable Headrest -----	104HE10	E0955	\$200
	<input type="checkbox"/> Links Multi Adjustable Curved Headrest -----	104HE11	E0955	\$240
	<input type="checkbox"/> Additional Link for Headrest -----	104HE12	E0955	\$50
	<input type="checkbox"/> 2 Additional Links for Headrest -----	104HE13	E0955	\$100
	<input type="checkbox"/> 3 Additional Links for Headrest -----	104HE14	E0955	\$150
	<input type="checkbox"/> UHRBA -----	104HEBA	E0955	\$61

FOOTNOTES:

- F. Kids Headrest Foam Pad available for chairs 10"W to 14"W
- G. Small Headrest Foam Pad available for chairs 15"W & 16"W
- H. Composite Footplates available for seat width greater than or equal to 14". For seat widths less than 14" use Super Kids Aluminum Flip Up footplates

19 HEADREST FOAM PAD OPTIONS

	<input type="checkbox"/> Kids Headrest Foam Pad-----	104HF5		Standard
	<input type="checkbox"/> Small Headrest Foam Pad -----	104HF1		No Charge
	<input type="checkbox"/> Small Headrest with Wings Foam Pad -----	104HF2		\$90
	<input type="checkbox"/> Medium Headrest Foam Pad -----	104HF7		No Charge
	<input type="checkbox"/> Medium Headrest with Wings Foam Pad -----	104HF8		\$90
	<input type="checkbox"/> Large Headrest Foam Pad -----	104HF3		No Charge
	<input type="checkbox"/> Large Headrest with Wings Foam Pad -----	104HF4		\$95
	<input type="checkbox"/> Omit Headrest -----	104HE09		less net \$10

20 ANTI TIP


	<input type="checkbox"/> Fixed Anti-Tipper -----	104MA22		\$75
--	--	---------	--	------

21 MISCELLANEOUS ACCESSORIES


	<input type="checkbox"/> Auto-Seat Positioning Straps-----	104MA03	E0971	\$52
	<input type="checkbox"/> Velcro Release Positioning Strap-----	104MA05	E0978	\$48
	<input type="checkbox"/> Quick Release Axles-----	104AL1		\$60
	<input type="checkbox"/> Padded Lap Tray-----	104MA09	E0950	\$285
	<input type="checkbox"/> Half Padded Lap Tray, Left-----	104MA14	E0950	\$160
	<input type="checkbox"/> Half Padded Lap Tray, Right-----	104MA15	E0950	\$160

SEAT TO FLOOR HEIGHT


STF FRONT (inches)	REAR WHEEL (inches)	CASTER SIZE (inches)	FORK SIZE	FORK HOLE	AXLE PLATE HOLE	STF REAR (inches)
18.50	12	4	LG	C	A	17.25
18.50	12	5	LG	B	A	17.25
19.00	12	7	SML	C	A	17.50
19.50	12	7	LG	B	B	18.25
20.00	12	8	LG	B	B	18.50
17.50	20	4	SML	C	A	15.75
18.00	20	5	SML	C	A	15.75
18.50	22	5	LG	B	A	16.75
19.00	20	7	SML	C	B	17.00
19.50	22	7	LG	B	B	18.00
20.00	22	7	LG	C	B	18.00




LONG FORK



SHORT FORK



**12" WHEEL
AXLE PLATE**



**20", 22" WHEEL
AXLE PLATE**